

Crime-lite AUTO® ACCESSORIES



Clip-on ring light attachment for improved imaging of tool marks and 3D fingerprints

Crime-lite AUTO

OBLIQUE Ring Light

PART# QCL/AUTO/ACC/9

A 24x LED clip-on ring light attachment for the Crime-lite AUTO that will improve ridge detail when capturing 3D fingermarks in dust, mud, dried blood etc.

Offering directional, high-intensity illumination via 12x White and 12x UV LEDs (software selectable), the Oblique Ring Light attachment provides multi-angle surface lighting at all angles around an item of evidence.

Ordering Information

CRIME-LITE AUTO - OBLIQUE RING LIGHT
QCL/AUTO/ACC/9

- 5700K White or UV 365nm (Software selectable)
- 12 LEDs per wavelength
- Intensity control
- Lighting direction control – Quadrant selection
- Powered via USB connection to Crime-lite AUTO



foster+freeman
Forensic Science Innovation

Head Office, UK Sales Office
Vale Park | Evesham
WR11 1TD | United Kingdom

T: +44 (0)1386 768 050
E: sales@fosterfreeman.com

USA Sales Office
46030 Manekin Plaza | Suite 170
Sterling | VA 20166 | USA

T: 888 445 5048
E: usoffice@fosterfreeman.com

EU Sales Office
Foster+Freeman GmbH
Düsseldorf | Germany

T: +49 (0)211 158 422 57
E: owen.lang@fosterfreeman.com

NL Sales Office
Foster+Freeman B.V.
Heerhugowaard | Nederland

T: +31 (0)6 1114 3858
E: nl.sales@fosterfreeman.com

Crime-lite AUTO[®] ACCESSORIES


Crime-lite AUTO

OBLIQUE Ring Light


PART# QCL/AUTO/ACC/9



Left LEDs Active




Plastic fingerprint in White Silicone




Right LEDs Active




Top and Left LEDs Active



Plastic fingerprint in Red Silicone



Right and Bottom LEDs Active



Left LEDs Active



Right LEDs Active




Left and Right LEDs Active


Toolmarks in Wood



Bottom LEDs Active



Top LEDs Active



Top and Bottom LEDs Active

foster+freeman
Forensic Science Innovation

Head Office, UK Sales Office
Vale Park | Evesham
WR11 1TD | United Kingdom

T: +44 (0)1386 768 050
E: sales@fosterfreeman.com

USA Sales Office
46030 Manekin Plaza | Suite 170
Sterling | VA 20166 | USA

T: 888 445 5048
E: usoffice@fosterfreeman.com

EU Sales Office
Foster+Freeman GmbH
Düsseldorf | Germany

T: +49 (0)211 158 422 57
E: owen.lang@fosterfreeman.com

NL Sales Office
Foster+Freeman B.V.
Heerhugowaard | Nederland

T: +31 (0)6 1114 3858
E: nl.sales@fosterfreeman.com